



Blue (/52) TL/TL-D

TL 20W/52 SLV/25

By emitting light almost entirely within the 400 to 500 nm bandwidth these Medical Therapy Jaundice TL/TL-D /52 lamps have no radiation from the short wave UVB waveband. They are therefore ideal for treating new born babies suffering from hyperbilirubinemia (neonatal jaundice) and Crigler-Najjar Syndrome (CNS). Moreover, the bandwidth of these lamps peak at the most effective treatment wavelength of 450 nm. This highly efficacious phototherapy treatment has eliminated the need for blood transfusions in almost all jaundiced infants.

Warnings and Safety

- A lamp breaking is extremely unlikely to have any impact on your health. If a lamp breaks, ventilate the room for 30 minutes and remove the parts, preferably with gloves. Put them in a sealed plastic bag and take it to your local waste facilities for recycling. Do not use a vacuum cleaner.

Product data


General Information		Voltage (Nom)	
Cap-Base	G13 [Medium Bi-Pin Fluorescent]	59 V	
Life to 50% Failures (Nom)	9,000 hour(s)	Controls and Dimming	
Useful Life (Nom)	2,000 hour(s)	Dimmable	No
Light Technical		Mechanical and Housing	
Color Code	52	Bulb Shape	T38
Luminous Flux	318 lm	Net Weight (Piece)	156.000 g
Chromaticity Coordinate X (Nom)	154	Approval and Application	
Chromaticity Coordinate Y (Nom)	69	Mercury (Hg) Content (Nom)	8.0 mg
Operating and Electrical		CE mark	Yes
Power Consumption	19.3 W	Product Data	
Lamp Current (Nom)	0.36 A	Order product name	TL 20W/52 SLV/25

Blue (/52) TL/TL-D

Full product name	TL 20W/52 SLV/25
Full product code	871150064302540
Order code	64302540
Material Nr. (12NC)	928003505203
Numerator - Quantity Per Pack	1


EAN/UPC - Product/Case	8711500643025
Numerator - Packs per outer box	25
EAN/UPC - Case	8711500643032

Dimensional drawing



Product	D (max)	A (max)	B (max)	B (min)	C (max)
TL 20W/52 SLV/25	40.5 mm	589.8 mm	596.9 mm	594.5 mm	604 mm

Photometric data



Spectral Power Distribution Colour - TL 20W/52 SLV/25

